

<b>Faecal Occult Bottles with and without dosing filter</b>
---

## Faecal Occult Bottles

For the manufacturers of test kits, for the diagnose of blood or stool, we have an extended range of Faecal Occult Bottles.

Composition: bottle and screw cap with mounted spatula and with “dropping-end”

- made from medical PP and LDPE
- Leak-proof design
- Different colour of caps



Cat. No	Volume	Height	Diameter	Total Height	Material	Colour	Packing
FOB-07	5 ml	47 mm	17 mm	65 mm	LDPE, PP	Natural + red	5000 pcs/ctn
FOB-08	5 ml	56 mm	15 mm	76 mm	LDPE, PP	Natural + blue	5000 pcs/ctn
FOB-09	5 ml	56 mm	13 mm	76 mm	LDPE, PP	Natural + purple	5000 pcs/ctn

## Faecal Occult Bottles with dosing-filter

- With dosing filter in order to take the same quantity of stool
- Made from medical grade LDPE and PP
- Screw cap with spatula



Cat. No	Volume	Height	Diameter	Total Height	Material	Colour	Packing
FOB-209	3 ml	88 mm	12 mm	106 mm	LDPE, PP	Natural + white/green	5000 pcs/ctn
FOB-205	5 ml	85 mm	13 mm	104 mm	LDPE, PP	Natural + red/blue	5000 pcs/ctn