




Reusable Glass Slides for Latex tests

- Use for latex and other agglutination tests
- Large wells
- Customized layout is possible
- Teflon printed ink
- 2 mm thick glass plate with ground edges (for optimal safety)
- Black, white or transparent rear side

Standard Glass Slides

Cat. No	Description	Size	Example
61.802.01/B	3 blue circular numbered rings (diam. 30 mm) Black rear side	52 x 108 mm	
61.802.02/B	6 BLUE oval rings (30 x 19 mm) Black rear side	52 x 108 mm	
61.802.02/R	6 RED oval rings (30 x 19 mm) Black rear side	52 x 108 mm	
61.802.02/G	6 GREEN oval rings (30 x 19 mm) Black rear side	52 x 108 mm	
61.802.02/W	6 WHITE oval rings (30 x 19 mm) Black rear side	52 x 108 mm	

61.801.02/B	6 BLUE oval rings (30 x 19 mm) White rear side	52 x 108 mm	
61.801.02/BL	6 BLACK oval rings (30 x 19 mm) White rear side	52 x 108 mm	

Examples customized Glass Slides

