



Pipettes

For use in diagnostic test kits we have a large selection of very accurate disposable plastic μ -pipettes for dispensing blood, serum, plasma, urine, or other fluids such as Latex.

The various drop sizes, bulb sizes and shapes of these pipettes are able to meet the different requirements for Diagnostic Tests such as Latex tests, RPR tests and Rapid Diagnostic Tests.

Several shapes and drop sizes

- The pipettes without a bulb are supposed to dispense only one single drop.
- Pipettes with a bulb can dispense several drops.
- Some pipettes have a paddle. With this paddle you can spread the liquid on the test card. This provides the advantage of dual-usage, as dispensing pipette and as a mixing stick.
- Pipettes are available with different drop size.
- Our standard range of pipettes will meet most of the requirements. But we can also customize the design following your requirements.

Customized packaging

In order to meet your needs we can customize the packaging. For your comfort we can pre-count these pipettes following your specifications and we supply them pre-packed in plastic bags with a zip-closure or in boxes or just with an elastic band.

Standard pipettes

We offer 4 types of disposable plastic μ -pipettes for use in diagnostic test kits

1) the "mixing" pipette :

- these pipettes dispense only one single drop
- the paddle end of the pipette can be used as mixing stick

2) the "bulb" pipette :

- these pipettes can dispense several drops

3) the "mixing-bulb" pipette :

- these pipettes can dispense several drops
- the paddle end of the pipette can be used as mixing stick

Cat. No.	Type	Drop size	Number of drops (+/-)	Length	Standard pre-packaging
P-00.50	Mixing pipette	50 μ l	1 drop	115 mm	
P-00.50/100	Mixing pipette	50 μ l	1 drop	115 mm	100/bag
P-04.40/VR	Bulb-pipette	40 μ l	15 drops	85 mm	2000/bag – 20000/box
P-05.20/S10Z	Bulb pipette	20 μ l	30 drops	105 mm	10/bag - sterile
P-10.30	Mixing-bulb pipette	30 μ l	6 drops	130 mm	500/box – 10000/box
P-10.50/100	Mixing bulb pipette	50 μ l	4 drops	125 mm	100/bag
P-10.50	Mixing bulb pipette	50 μ l	4 drops	125 mm	