

## Possibilities

- Different material
  - o Laminated cardboard with mat or shiny surface (black or white)
  - o Mat PVC card
  - o Shiny PVC card
- Can be embossed
- We can make “tear-off” lines for splitting the card easily
- Plasma H treatment of the card
- Hydrophobic surface with no cross contamination between the wells
- Perfect spreading properties
- Different printing processes
  - o Conventional Screen Printing
  - o UV Screen Printing
  - o Offset Printing
- Packaging
  - o Bulk : cards are packed per 100 cards. On top a transparent PVC protection card, on the bottom a grey cardboard card and packed by an elastic band
  - o Pre-count: counted cards in packs, packed by an elastic band
  - o Pre-Packed: Counted cards packed in a Cellophane bag or Zip-lock bag
- Standard and Customized cards

**Please ask our help for receiving samples and all possibilities for your customized card!**

**Please find our standard cards on the next page**









Pre- counted cards



Packed cards

## Standard Cards

Cat. No	Description	Card-size	Material	Example	Packing
61.815.00	6 black mat printed circular fields	60 x 95 mm	White shiny PVC		5000 pcs/ctn
61.816.00	6 white mat printed circular fields with red ring	60 x 95 mm	White shiny PVC		5000 pcs/ctn
61.890.10	RPR card – 10 embossed circles (18 mm) with red ring	75 x 124 mm	Laminated cardboard with white mat film		5000 pcs/ctn
61.898.10	RPR card – not embossed – 10 grey printed	75 x 124 mm	Laminated cardboard with white mat film		5000 pcs/ctn

	circles (18 mm)				
61.821.02	6 green oval rings with white 'write-on' fields	75 x 118 mm	Laminated cardboard with black mat film		5000 pcs/ctn
61.820.02	6 red oval rings and white 'writ-on' fields	75 x 118 mm	Laminated cardboard with black shiny film		5000 pcs/ctn

## Customized cards

For a production (or quotation) for custom-made cards we need following information:

- Material of the card
  - o White PVC
  - o Black PVC
  - o Laminated cardboard with a white mat film
  - o Laminated cardboard with a black mat film
  - o Laminated cardboard with a black shiny film
  - o ...
- Colour of the surface of the card
- Printed circles or printed fields
- Design
- Card size
- Packaging
  - o Bulk
  - o Pre-count
  - o Packed in cellophane packaging or zip-lock bag
- Quantity per production
  - o Depends of the type of card, in general minimum order of customized cards between 20.000 – 50.000 cards

