

Type PSZ

- 4 ml up to 100 ml vials
- In natural or Black colour
- Different colour of caps for colour coding: yellow, red, blue, white, green, black
- Made from medical PP and HDPE
- Leak-free closure
- Resistant to regular chemical solution
- Bottles 4 ml to 6 ml have 15 mm screw neck
- Bottles 10 ml to 100 ml have 22 mm screw neck



Bottles

Cat. No	Volume	Height	Diameter	Material	Colour	Packing
PSZ-004N	4 ml	32 mm	18 mm	HDPE	Natural	5000 pcs/ctn
PSZ-004BL	4 ml	32 mm	18 mm	HDPE	Black	5000 pcs/ctn
PSZ-006N	6 ml	39 mm	18 mm	HDPE	Natural	3000 pcs/ctn
PSZ-006BL	6 ml	39 mm	18 mm	HDPE	Black	3000 pcs/ctn
PSZ-010N	10 ml	48 mm	25 mm	HDPE	Natural	1000 pcs/ctn
PSZ-010BL	10 ml	48 mm	25 mm	HDPE	Black	1000 pcs/ctn
PSZ-020N	20 ml	60 mm	27 mm	HDPE	Natural	1000 pcs/ctn
PSZ-020BL	20 ml	60 mm	27 mm	HDPE	Black	1000 pcs/ctn
PSZ-030N	30 ml	60 mm	32 mm	HDPE	Natural	1000 pcs/ctn
PSZ-050N	50 ml	61 mm	43 mm	HDPE	Natural	1000 pcs/ctn
PSZ-050BL	50 ml	61 mm	43 mm	HDPE	Black	1000 pcs/ctn
PSZ-100N	100 ml	102 mm	43 mm	HDPE	Natural	500 pcs/ctn

Caps

Cat. No	Colour	Material	Information	Packing
PSZ-C15-Y	Yellow	PP	For PSZ-bottles 4 to 6 ml	10000 pcs/ctn
PSZ-C15-B	Blue	PP	For PSZ-bottles 4 to 6 ml	10000 pcs/ctn
PSZ-C15-BL	Black	PP	For PSZ-bottles 4 to 6 ml	10000 pcs/ctn
PSZ-C15-O	Orange	PP	For PSZ-bottles 4 to 6 ml	10000 pcs/ctn
PSZ-C15-G	Green	PP	For PSZ-bottles 4 to 6 ml	10000 pcs/ctn
PSZ-C15-W	White	PP	For PSZ-bottles 4 to 6 ml	10000 pcs/ctn
PSZ-C15-R	Red	PP	For PSZ-bottles 4 to 6 ml	10000 pcs/ctn
PSZ-C22-Y	Yellow	PP	For PSZ-bottles 10 to 100 ml	10000 pcs/ctn
PSZ-C22-B	Blue	PP	For PSZ-bottles 10 to 100 ml	10000 pcs/ctn
PSZ-C22-BL	Black	PP	For PSZ-bottles 10 to 100 ml	10000 pcs/ctn
PSZ-C22-O	Orange	PP	For PSZ-bottles 10 to 100 ml	10000 pcs/ctn
PSZ-C22-G	Green	PP	For PSZ-bottles 10 to 100 ml	10000 pcs/ctn
PSZ-C22-W	White	PP	For PSZ-bottles 10 to 100 ml	10000 pcs/ctn
PSZ-C22-R	Red	PP	For PSZ-bottles 10 to 100 ml	10000 pcs/ctn

Type PSN



Reagent bottles

Cat. No	Volume	Height	Diameter	Material	Colour	Packing
PSN-015N/x	15 ml	58 mm	25 mm	HDPE	Natural	1000 pcs/ctn
PSN-030N/x	30 ml	61 mm	34 mm	HDPE	Natural	1000 pcs/ctn
PSN-060N/x	60 ml	84 mm	38 mm	HDPE	Natural	1000 pcs/ctn
PSN-015A/x	15 ml	58 mm	25 mm	HDPE	Amber	1000 pcs/ctn
PSN-030A/x	30 ml	61 mm	34 mm	HDPE	Amber	1000 pcs/ctn
PSN-060A/x	60 ml	84 mm	38 mm	HDPE	Amber	1000 pcs/ctn