

Type PSZ

- 4 ml up to 100 ml vials
- In natural or Black colour
- Different colour of caps for colour coding: yellow, red, blue, white, green, black
- Made from medical PP and HDPE
- Leak-free closure
- Resistant to regular chemical solution
- Bottles 4 ml to 6 ml have 15 mm screw neck
- Bottles 10 ml to 100 ml have 22 mm screw neck



Bottles

| Cat. No | Volume | Height | Diameter | Material | Colour | Packing |
|-----------|--------|--------|----------|----------|---------|--------------|
| PSZ-004N | 4 ml | 32 mm | 18 mm | HDPE | Natural | 5000 pcs/ctn |
| PSZ-004BL | 4 ml | 32 mm | 18 mm | HDPE | Black | 5000 pcs/ctn |
| PSZ-006N | 6 ml | 39 mm | 18 mm | HDPE | Natural | 3000 pcs/ctn |
| PSZ-006BL | 6 ml | 39 mm | 18 mm | HDPE | Black | 3000 pcs/ctn |
| PSZ-010N | 10 ml | 48 mm | 25 mm | HDPE | Natural | 1000 pcs/ctn |
| PSZ-010BL | 10 ml | 48 mm | 25 mm | HDPE | Black | 1000 pcs/ctn |
| PSZ-020N | 20 ml | 60 mm | 27 mm | HDPE | Natural | 1000 pcs/ctn |
| PSZ-020BL | 20 ml | 60 mm | 27 mm | HDPE | Black | 1000 pcs/ctn |
| PSZ-030N | 30 ml | 60 mm | 32 mm | HDPE | Natural | 1000 pcs/ctn |
| PSZ-050N | 50 ml | 61 mm | 43 mm | HDPE | Natural | 1000 pcs/ctn |
| PSZ-050BL | 50 ml | 61 mm | 43 mm | HDPE | Black | 1000 pcs/ctn |
| PSZ-100N | 100 ml | 102 mm | 43 mm | HDPE | Natural | 500 pcs/ctn |

Caps

| Cat. No | Colour | Material | Information | Packing |
|------------|--------|----------|------------------------------|---------------|
| PSZ-C15-Y | Yellow | PP | For PSZ-bottles 4 to 6 ml | 10000 pcs/ctn |
| PSZ-C15-B | Blue | PP | For PSZ-bottles 4 to 6 ml | 10000 pcs/ctn |
| PSZ-C15-BL | Black | PP | For PSZ-bottles 4 to 6 ml | 10000 pcs/ctn |
| PSZ-C15-O | Orange | PP | For PSZ-bottles 4 to 6 ml | 10000 pcs/ctn |
| PSZ-C15-G | Green | PP | For PSZ-bottles 4 to 6 ml | 10000 pcs/ctn |
| PSZ-C15-W | White | PP | For PSZ-bottles 4 to 6 ml | 10000 pcs/ctn |
| PSZ-C15-R | Red | PP | For PSZ-bottles 4 to 6 ml | 10000 pcs/ctn |
| PSZ-C22-Y | Yellow | PP | For PSZ-bottles 10 to 100 ml | 10000 pcs/ctn |
| PSZ-C22-B | Blue | PP | For PSZ-bottles 10 to 100 ml | 10000 pcs/ctn |
| PSZ-C22-BL | Black | PP | For PSZ-bottles 10 to 100 ml | 10000 pcs/ctn |
| PSZ-C22-O | Orange | PP | For PSZ-bottles 10 to 100 ml | 10000 pcs/ctn |
| PSZ-C22-G | Green | PP | For PSZ-bottles 10 to 100 ml | 10000 pcs/ctn |
| PSZ-C22-W | White | PP | For PSZ-bottles 10 to 100 ml | 10000 pcs/ctn |
| PSZ-C22-R | Red | PP | For PSZ-bottles 10 to 100 ml | 10000 pcs/ctn |

Type PSN



Reagent bottles

| Cat. No | Volume | Height | Diameter | Material | Colour | Packing |
|------------|--------|--------|----------|----------|---------|--------------|
| PSN-015N/x | 15 ml | 58 mm | 25 mm | HDPE | Natural | 1000 pcs/ctn |
| PSN-030N/x | 30 ml | 61 mm | 34 mm | HDPE | Natural | 1000 pcs/ctn |
| PSN-060N/x | 60 ml | 84 mm | 38 mm | HDPE | Natural | 1000 pcs/ctn |
| PSN-015A/x | 15 ml | 58 mm | 25 mm | HDPE | Amber | 1000 pcs/ctn |
| PSN-030A/x | 30 ml | 61 mm | 34 mm | HDPE | Amber | 1000 pcs/ctn |
| PSN-060A/x | 60 ml | 84 mm | 38 mm | HDPE | Amber | 1000 pcs/ctn |